


The Ascent Camera Series supports a variety of interline transfer sensors with high quantum efficiency and low dark current. Most cameras are available with a choice of monochrome sensor or color sensor with a Bayer mosaic. Low noise and small pixels are ideal for chemiluminescence, fluorescence, and microscopy.

- CCDs up to 29 megapixels
- 16-bit digitization up to 16 MHz throughput
- Dual channel CCD support
- Video focus mode
- USB 2.0 interface: no plug in cards or external controllers
- Programmable, intelligent cooling to 35°C below ambient
- Binning up to 8 Horizontal x CCD height
- Subarray readout and fast sequencing modes
- Programmable offset and gain
- External triggering and strobe controls
- ActiveX drivers included with every system
- Field upgradeable firmware
- Fused silica window
- Optional 6-position 1" or 8-position 1.25" filter wheels
- Optional C-mount, Nikon F-mount, or 2" slip fit adapter
- Single 6V supply
- Compact enclosure: 23 oz. (0.65 kg)
- Programmable status indicators



Clockwise from lower left: image courtesy Scientific Instrument Company; slide courtesy Ventana Medical Systems; image courtesy Dr. David Rapaport, UCSD

CCD SPECIFICATIONS

Model	CCD	Array	Pixels	Pixel size (microns)	CCD Size (mm)	Area (mm ²)	Diagonal (mm)	Video size (")
A29050	KAI-29050	6576 x 4384	28.8M	5.5	36 x 24	872	43.5	2.7
A16000	KAI-16000	4872 x 3248	15.8M	7.4	36 x 24	866.5	43.3	2.7
A8050	KAI-8050	3296 x 2472	8.1M	5.5	18.1 x 13.6	246.5	22.7	1.42
A4000	KAI-4022	2048 x 2048	4.2M	7.4	15.2 x 15.2	229.7	21.4	1.34
A4050	KAI-4050	2336 x 1752	4.2M	5.5	11.3 x 11.3	126.9	15.9	1.0
A2150	KAI-2150	1920 x 1080	2.1M	5.5	10.6 x 5.9	62.7	12.1	0.76
A2000	KAI-2020	1600 x 1200	1.9M	7.4	11.8 x 8.9	105.1	14.8	0.93
A2050	KAI-2050	1600 x 1200	1.9M	5.5	8.8 x 6.6	58.1	11.0	0.69
A1050	KAI-1050	1024 x 1024	1.1M	5.5	5.6 x 5.6	34.8	8.3	0.52
A285	ICX285	1360 x 1024	1.4M	6.45	8.8 x 6.6	57.9	11.0	0.69
A205	ICX205	1360 x 1024	1.4M	4.65	6.3 x 4.8	30.1	7.9	0.49
A340	KAI-0340	684 x 484	313K	7.4	4.8 x 3.6	17.2	6.0	0.37

For complete CCD specifications, including cosmetic grading, see data sheet from manufacturer.

CAMERA DATA SHEETS

A complete set of camera data sheets as well as mechanical drawings are on our Integration Starter Kit CD, or at www.ccd.com



151 N. Sunrise Ave. Ste 902
Roseville CA 95661 USA
tel 916 218 7450
fax 916 218 7451
www.ccd.com



PC Interface	USB 2.0
USB2 Cable	Std.: 5m. Extensions: 5 meters between hubs; 5 hubs maximum (max. total of 30m) Wide variety of extenders available, including fiber optics to 10 km.
Digital Resolution	16 bits at up to 16 MHz throughput*
System Noise (example)	4.5 e ⁻ RMS at 3.3 MHz (A285); 6 e ⁻ RMS at 10 MHz (A285); 7 e ⁻ RMS at 8/16 MHz (A2150)
Pixel Binning	1 x 1 to 8 x height of CCD on-chip
Exposure Time	Minimum 100 microseconds; max.183 minutes
Image Sequencing	1 to 65535 image sequences under software control
Frame Sizes	Full frame, subframe, focus mode
Cooling (typical)	Thermoelectric cooler. Maximum forced air cooling to 40°C below ambient temperature.
Temperature Stability	± 0.1°C
Camera Head Size	Aluminum. 3.2" x 4.7" x 1.3" (8.1 x 11.9 x 3.3 cm) Weight: 1.4 lb. (0.65 kg)
Mounting	1.5" x 2.5" bolt pattern, 6-32 thread. Optional C-mount (1" 32 tpi thread), Nikon F-mount adapters.
Back Focal Distance	Standard: 0.32" (0.81 cm). [optical]
Operating Environment	-30°C to 35°C. Relative humidity: 10 to 90% non-condensing.
Op.Sys.Support	Windows, Linux, Mac OSX
Power	20W maximum power with internal shutter open and cooling maximum. AC/DC "brick" supply with int'l AC input plug (100-240V, 50-60 Hz). Alternate 6V input from user's source.
Remote Triggering	LVTTL input allows exposure to start within 25 microseconds of rising edge of trigger

*Varies from CCD to CCD. See data sheet for specific model.

CCD SENSITIVITY

